

Using FT-IR Spectroscopy to Elucidate the Structures of Ablative Polymers

Wendy Fan

The composition and structure of an ablative polymer has a multifaceted influence on its thermal, mechanical and ablative properties. Understanding the molecular level information is critical to the optimization of material performance because it helps to establish correlations with the macroscopic properties of the material, the so-called structure-property relationship. Moreover, accurate information of molecular structures is also essential to predict the thermal decomposition pathways as well as to identify decomposition species that are fundamentally important to modeling work. In this presentation, I will describe the use of infrared transmission spectroscopy (FT-IR) as a convenient tool to aid the discovery and development of thermal protection system materials.

Using FT-IR Spectroscopy to Elucidate the Structures of Ablative Polymers

Wendy Fan

ERC Inc.

NASA Ames Research Center



Outline

1. Introduction

- Importance of polymer structure (Motivation)
- FT-IR principle and technique


2. Applications

- Detailed structural analysis of a phenolic polymer
- Correlation among cure condition - crosslinking - thermal stability/char yield
- Correlation between crosslinking - thermal oxidative decomposition mechanism
- Diagnostic of a thermal decomposition product

3. Conclusion



Structure of Hemoglobin – How Nature Designs Its Materials



Unique
sequence of
amino acids

Held together
by hydrogen
bonding

Held together
by hydrogen,
ionic, sulfide
bonding


Several
peptide chains
joined together

<http://alevelnotes.com>



Polymer Structure and Morphology


...influence key thermal, mechanical and ablative properties



Good morphology



Poor morphology





Techniques to Identify Polymer Structures

1. Nuclear Magnetic Resonance (NMR, both liquid and solid state)

- Precise information of how atoms are bonded

2. FT-IR (Fourier Transform Infrared) transmission spectroscopy

- Identify polymer structures, in particular polar functional groups, crosslink degree, crystalline region

3. Raman spectroscopy


- Identify non-polar groups, such as olefinic groups (as in graphitic species)

4. Mass spectrometry


- Information about fragments/groups, fragmentation modes

5. XPS (x-ray Photoelectron Spectroscopy)


- Information about elements and functional groups on a surface



FT-IR – An Introduction



<http://www.biophysik.uni-freiburg.de/Spectroscopy/FTIR/spectroscopy.html>



Main advantages of FT-IR

- Speed
- Sensitivity
- High spectral accuracy and resolution
- Internally calibrated
- Mechanical simplicity
- Easy sample preparation and operation




andersonmaterials.com




FT-IR spectroscopy: Theory


Vibrational modes of a -CH₂ group




Symmetrical stretching




Twisting




Wagging



Scissoring



Rocking



Asymmetric stretching

Vibrational transitions absorb discrete energy levels unique to a molecule


http://en.wikipedia.org/wiki/Fourier_transform_infrared_spectroscopy



FT-IR Spectrum of Formaldehyde



C




- To chemists, a FT-IR spectrum is like a pair of eyes looking into a molecule
- The intensity, shape, width and location of a peak all reflect the specific local environment the group is in



Curing Chemistry of Phenolic Resin


Phenolic resin (resole)




Cure



Phenolic Network polymer




A typical reaction in network formation







FT-IR of a Fully Crosslinked Phenolic Polymer



- Textbook example and extremely informative
- Absorption peaks can be easily assigned to corresponding groups
- Quantitative estimation of the degree of crosslinking is possible.





Region 1: 4000cm^{-1} - 2700cm^{-1}



3500 cm^{-1} : Phenolic OH
3040 cm^{-1} : aromatic sp^2 C-H
2940 cm^{-1} : Aliphatic sp^3 CH_2 asym. Stretch
2860 cm^{-1} : Aliphatic sp^3 CH_2 sym. stretch



Region 2: 1700cm^{-1} - 1000cm^{-1}



$1610\text{-}1650\text{cm}^{-1}$: Benzene $\text{C}=\text{C}$

1475cm^{-1} : CH_2 scissor

1260cm^{-1} : CH_2 rocking

1210cm^{-1} : sp^2 $\text{C}-\text{O}$







1160cm^{-1} : CH_2 wagging

1120cm^{-1} : sp^2 $\text{C}-\text{H}$

$I_{\text{CH}_2} / I_{\text{C}=\text{C}}$ can be used to semi-quantitatively measure the crosslinking degree



Region 3: 900cm⁻¹ - 700cm⁻¹




Most highly crosslinked

13
1,2,4,6-ortho, para, ortho-



Application




Experiments: phenolic resin was cured under three different conditions, each with increasing temperature and cure time.

Condition 1, Condition 2, Condition 3

→
Temperature and cure time increases




FT-IR of Cured Phenolic Polymers



Change of intensity in multiple peaks indicates the progress of curing (more crosslink




Thermogravimetric Analysis (TGA) (in Nitrogen, 20°C/min.)



Char yield correlates with degree of polymer crosslinking




TGA in Air 20°C/min.



Weight% increases with increased crosslinking



Decomposition of Phenolic – Char Forming Mechanism




Char structures were
characterized by solid state ^{13}C
NMR


Substructures of graphitic/carbon char

**Virgin polymer structure dictates the
decomposition pathways and char yield**

Amorphous carbon char



Application 2: Identifying a Decomposition Process




- $-\text{CH}_3$ group formed
- Lower degree of crosslinking



Thermal Decomposition in Vacuum

Likely Pathway




- Lower degree of crosslinking
- -CH_3 group formed

G. E. Maciel and I.-S. Chuang, *Macromolecules*, 1984, 17, 1081-1087
Fyfe, C. A. et al, *Macromolecules*, 1983, 16, 1216-1219
Zinke, A. *J. Appl. Chem.*, 1951, 1, 257.



FT-IR Spectrum of a Model Polymer



Center portion


Surface portion

Model polymer

Model polymer was cured in nitrogen and at an elevated temperature (close to the oven condition)



Bottom Up Approach to TPS Materials Discovery



- A thorough understanding of chemistry/molecular structure is key
- Most of these correlations have been obtained by experiments



Advantages and Limitations

Advantages:


1. Convenient, low cost, speedy
2. Excellent tool for first screening of polymer structures and identification of unknowns
3. Useful for predicting reactivity and physicochemical properties, as well as interpreting decomposition mechanisms and products
4. Routinely used at Ames for qualitative determination of new polymer structures as well as analyzing curing products
5. Can be coupled with TGA to study thermal decomposition products (TGA-IR)

Limitations:

1. Mostly qualitative, does not provide precise molecular structures
2. Can be semi-quantitative with the right software
3. Not all polymers have well-resolved absorption peaks in the IR region (overlaps, weak absorptions)





Using TGA-FT IR to Monitor Gas Phase Products



http://www.uclan.ac.uk/schools/forensic_investigative/fire_hazards_science/equipment_tests/tga_ftir_magna.php

- TGA-IR allows the gas phase products to be collected continuously as a function of temperature.
- Gas phase infrared spectra are much sharper than condensed phase spectra, and recognition of the individual molecules is possible.



<http://cnx.org/content/m23038/latest/>

Contributors

Tane Boghozian, Ehson Ghandehari, Firouzeh
Mohadjerani, Jeremy Thornton

Acknowledgement

Hypersonics group of
Thermal Protection Materials branch,
NASA Ames Research Center

Funding

Hypersonics
EDL-TDP